

## Symbolerklärung

Auf den Etiketten der Dentaco Produkte werden verschiedene Piktogramme verwendet.



brennbar



Achtung



Gesundheitsgefahr



Dieses Produkt trägt das CE-Kennzeichen in Übereinstimmung mit den Bestimmungen der Richtlinie 93/42/EWG des Rates vom 14. Juni 1993 über Medizinprodukte.



Chargenbezeichnung



Informationsbeilage bzw. Hinweise des Herstellers beachten.

XXX ml Nennfüllmenge e maximal zulässige Abweichungen werden eingehalten



Randvollvolumen



Produkt ist nach Aerosolverpackungsordnung konform



maximale Temperatur bei Lagerung und Transport: 50°C



Verfalldatum



Artikelnummer



Observe the instructions for use

## Symbol explanation

Various pictograms are used on the labels of Dentaco products.



flammable



harmful



Health hazard



This product bears the CE mark in accordance with the provisions of Council Directive 93/42/EEC of 14 June 1993 on medical devices.



batch number



Observe accompanying documents and instructions of manufacturer

XXX ml nominal filling quantity



brim full volume



article number



according to aerosol packing regulation



maximum temperature during storage and transport: 50°C



maximum permissible deviations are observed



Expiry date



Observe the instructions for use



(Hersteller / manufacturer) DENTACO GmbH & Co. KG, Max-Keith-Str. 46; 45136 Essen (Germany)  
Phone: +49 (0) 201 / 802899-0; Fax: +49 (0) 201 / 802899-99; Mail: info@dentaco.de; Internet: www.dentaco.de